



A **unique** monoblock design

intra-operatively offering total **modularity**

- CCD 135°
- Offset 44mm
- Trunnion 12/14
- Narrow neck
- Material Titanium 6/4ELI
- Coatings – Bi-Coat (CP Ti. & HA)
- Distal Slots – Bi-directional (310mm stem only)

Proximal Screw



Distal Screw

Screws are available in Ø5mm & Ø7mm, lengths 20 -60mm in 5mm increments

Hip stems are available in 6 diameters ranging from Ø12mm - Ø22mm in 2mm increments in 4 lengths 210mm, 240mm, 270mm and 310mm



BAROS™ Revison Stem



Ø16 X 310mm
Stem




Ø16 X 210mm
Stem



Femoral Head		Offset	Leg Length	Lateral Offset
28mm	SN	-3.5mm	54.4	41.5
	MN	0	57.0	44.0
	LN	+3.5	59.5	46.5
	XLN	+7.0	61.8	48.8
	XXLM	+10.5	64.2	51.3
32mm	SN	-4.0	54.1	42.2
	MN	0	57.0	44.0
	LN	+4.0	59.8	46.8
	XLN	+8.0	62.6	49.7
	XXLN	+12.0	65.4	52.5

Proximal Screw & Washer



For Stems

- Ø5mm 14 & 16mm
- Ø7mm 18, 20 & 22mm

Distal Screw



For Stems

- Ø5mm 14 & 16mm
- Ø7mm 18, 20 & 22mm

March 10

DISTRIBUTED BY



Manufactured in the UK

Implants International Ltd
 71 Jay Avenue, Teesside Industrial Estate
 Thornaby-on-Tees, TS17 9LZ, United Kingdom
 tel: +44(0)1642 769080 fax: +44(0)1642 765848
 e-mail: enquiries@implantsinternational.com
 web: www.implantsinternational.com