




# IMPLANTS

INTERNATIONAL

ARTHROPLASTY



**CLST Uncemented**  
Rough Cermaic Blasted  
Hip Stem 9.00mm x 150mm

**CLST Unemented**  
CP Plasma Coated  
Hip Stem 7.00mm x 142.8mm


**CLST**  
**UNCEMENTED HIP STEM**



### CLST Hip Stem Uncemented CP Plasma

Stem Size (W)	Stem Length (L)	Taper	Cat#
7.00	142.8	12/14	HS1393
8.00	146.4	12/14	HS1394
9.00	150.0	12/14	HS1395
10.00	153.5	12/14	HS1396
11.25	158.0	12/14	HS1397
12.50	162.5	12/14	HS1398
15.00	171.5	12/14	HS1399
17.50	180.5	12/14	HS1400

Material: Titanium Alloy with CP Plasma Coating



### CLST Hip Stem Uncemented Rough Ceramic Blasted Finish

Stem Size (W)	Stem Length (L)	Taper	Cat#
7.00	142.8	12/14	HS1401
8.00	146.4	12/14	HS1402
9.00	150.0	12/14	HS1403
10.00	153.5	12/14	HS1404
11.25	158.0	12/14	HS1405
12.50	162.5	12/14	HS1406
15.00	171.5	12/14	HS1407
17.50	180.5	12/14	HS1408

Material: Titanium Alloy with Rough Ceramic Blasted Finish

April 10



Manufactured in the UK

DISTRIBUTED BY

Implants International Ltd  
 71 Jay Avenue, Teesside Industrial Estate  
 Thornaby-on-Tees, TS17 9LZ, United Kingdom  
 tel: +44(0)1642 769080 fax: +44(0)1642 765848  
 e-mail: enquiries@implantsinternational.com  
 web: www.implantsinternational.com