




IMPLANTS

INTERNATIONAL

ARTHROPLASTY



Muller Ribbed Straight
Hip Stem 15.0mm x 152mm

Muller Ribbed Straight
Hip Stem 6.35mm x 132mm


**MULLER RIBBED STRAIGHT
HIP STEM**



Muller Ribbed Straight Hip Stem

Stem Size (W)	Stem Length (L)	Taper	Cat#
5.00	127	12/14	HS1409
6.25	132	12/14	HS1410
7.50	137	12/14	HS1411
10.0	142	12/14	HS1412
12.5	147	12/14 </td <td>HS1413</td>	HS1413
15.0	152	12/14	HS1414

Material: Cobalt Chrome to ISO 5832-4



April 10

DISTRIBUTED BY



Manufactured in the UK

Implants International Ltd
 71 Jay Avenue, Teesside Industrial Estate
 Thornaby-on-Tees, TS17 9LZ, United Kingdom
 tel: +44(0)1642 769080 fax: +44(0)1642 765848
 e-mail: enquiries@implantsinternational.com
 web: www.implantsinternational.com