



# IMPLANTS

I N T E R N A T I O N A L

A R T H R O P L A S T Y

**Performa™** Femoral Heads  
in Wrought Cobalt Chrome



**Tribol™** Ceramic Heads  
in Aluminium Oxide



**OMNI-FIT™**  
FEMORAL HEADS



# OMNI-FIT™ FEMORAL HEADS

## Tribol™ Ceramic Heads in Aluminium Oxide

					Size	Offset	Cat #
Ø28mm					Ø28mm	SN	FH 3087
	SN -3.5	MN 0	LN +3.5		Ø28mm	MN	FH 3089
					Ø28mm	LN	FH 3088
Ø32mm					Ø32mm	SN	FH 3094
	SN -4.0	MN 0	LN +4.0	XLN +8.0	Ø32mm	MN	FH 3095
					Ø32mm	LN	FH 3096
					Ø32mm	XLN	FH 3097
Ø36mm					Ø36mm	SN	FH 3110
	SN -4.0	MN 0	LN +4.0	XLN +8.0	Ø36mm	MN	FH 3111
					Ø36mm	LN	FH 3112
					Ø36mm	XLN	FH 3113

## Performa™ Femoral Heads in Wrought Cobalt Chrome

Ø28mm						Ø28mm	SN	FH 3050
	SN -3.5	MN 0	LN +3.5	XLN +7.0	XXLN +10.5	Ø28mm	MN	FH 3051
						Ø28mm	LN	FH 3052
						Ø28mm	XLN	FH 3053
						Ø28mm	XXLN	FH 3054
Ø32mm						Ø32mm	SN	FH 3067
	SN -4.0	MN 0	LN +4.0	XLN +8.0	XXLN +12.0	Ø32mm	MN	FH 3068
						Ø32mm	LN	FH 3069
						Ø32mm	XLN	FH 3070
						Ø32mm	XXLN	FH 3104

March 10



Manufactured in the UK

DISTRIBUTED BY

Implants International Ltd  
 71 Jay Avenue, Teesside Industrial Estate  
 Thornaby-on-Tees, TS17 9LZ, United Kingdom  
 tel: +44(0)1642 769080 fax: +44(0)1642 765848  
 e-mail: enquiries@implantsinternational.com  
 web: www.implantsinternational.com