



# IMPLANTS

INTERNATIONAL

ARTHROPLASTY



## PERFORMA™

**Trapped Head Design:** 22.225mm  
**Taper Size:** 10/12mm  
**Cup Size:** 40mm - 54mm in 2mm increments  
**Neck Size:** Short, medium and long

**Trapped Head Design:** 26mm  
**Taper Size:** 12/14  
**Cup Size:** 42mm & 44mm  
**Neck Size:** Short, medium and extra long  
 XXXLN available as specials

**Trapped Head Design:** 28mm  
**Taper Size:** 12/14mm  
**Cup Size:** 46mm - 58mm in 2mm increments  
**Neck Size:** Short, medium, long and extra long  
 XXXLN available as specials

## PERFORMA S-FIT™

**Trapped Head Design:** 28mm  
**S-Fit™ Taper Size:** 11/13mm  
**Cup Size:** 46mm - 58mm in 2mm increments  
**Neck Size:** Short, medium, long and extra long  
 XXXLN available as specials

Detailed Catalogue numbers Available on request in electronic or hard copy format.

Material - Heads: Chrome Cobalt or Stainless Steel  
 Inserts: UHMWPE (Low Calcium Stearate)

PERFORMA™ BI-POLAR  
 &  
 PERFORMA S-FIT™ BI-POLAR



## 28MM Ø PERFORMA™ MODULAR SYSTEM



<b>Non-Trapped Head Design:</b>	28mm
<b>Taper Size:</b>	12/14mm
<b>Cup Size:</b>	46m – 58mm in 2mm increments
<b>Neck Size:</b>	Short, medium, long and extra long XXXLN available as specials

Detailed Catalogue numbers Available on request in electronic or hard copy format.

Material – Heads: Chrome Cobalt or Stainless Steel  
Inserts: UHMWPE (Low Calcium Stearate)

**INSTRUMENT KIT: 1 TRAY SYSTEM**  
**SURGICAL TECHNIQUE: REF ST# 016/99**

March 10

DISTRIBUTED BY



0120



Manufactured in the UK

Implants International Ltd  
71 Jay Avenue, Teesside Industrial Estate  
Thornaby-on-Tees, TS17 9LZ, United Kingdom  
tel: +44(0)1642 769080 fax: +44(0)1642 765848  
e-mail: [enquiries@implantsinternational.com](mailto:enquiries@implantsinternational.com)  
web: [www.implantsinternational.com](http://www.implantsinternational.com)