



IMPLANTS  
INTERNATIONAL  
ARTHROPLASTY



Titanium Acetabular Shell



Acetabular Liners

RIM-FIX™  
ACETABULAR SYSTEM



# RIM-FIX™ SYSTEM

## CP TITANIUM PLASMA COATED SHELL - 28MM

Size	Catalogue #
46mm	AS4001
48mm	AS4002
50mm	AS4003
52mm	AS4004
54mm	AS4005
56mm	AS4006
58mm	AS4007
60mm	AS4008
62mm	AS4009
64mm	AS4010
66mm	AS4011

MATERIAL: TITANIUM ALLOY- Ti coated

## STD LINER

Size	Head Size		
	28mm	32mm	36mm
46mm	AL3869	AL3870	-
48mm	AL3750	AL3758	-
50mm	AL3751	AL3759	AL3873
52mm	AL3752	AL3760	AL3874
54mm	AL3753	AL3761	AL3875
56mm	AL3754	AL3762	AL3876
58mm	AL3755	AL3763	AL3877
60mm	AL3756	AL3764	AL3878
62mm	AL3757	AL3765	AL3879
64mm	AL3861	AL3863	AL3880
66mm	AL3862	AL3864	AL3881

MATERIAL: UHMWPE

Catalogue #

## CP TITANIUM / HA PLASMA COATED SHELL - 28MM

Size	Catalogue #
46mm	AS4012
48mm	AS4013
50mm	AS4014
52mm	AS4015
54mm	AS4016
56mm	AS4017
58mm	AS4018
60mm	AS4019
62mm	AS4020
64mm	AS4021
66mm	AS4022

MATERIAL: TITANIUM ALLOY- Ti & HA coated

## LINER - 10 DEGREE HOOD

Size	Head Size		
	28mm	32mm	36mm
46mm	AL3871	AL3872	-
48mm	AL3766	AL3774	-
50mm	AL3767	AL3775	AL3882
52mm	AL3768	AL3776	AL3883
54mm	AL3769	AL3777	AL3884
56mm	AL3770	AL3778	AL3885
58mm	AL3771	AL3779	AL3886
60mm	AL3772	AL3780	AL3887
62mm	AL3773	AL3781	AL3888
64mm	AL3865	AL3867	AL3889
66mm	AL3866	AL3868	AL3890

MATERIAL: UHMWPE

Catalogue #

Instrument Kit: 1 Tray System

Surgical Technique: Ref ST# 015/99

April 10



Manufactured in the UK

DISTRIBUTED BY

Implants International Ltd  
 71 Jay Avenue, Teesside Industrial Estate  
 Thornaby-on-Tees, TS17 9LZ, United Kingdom  
 tel: +44(0)1642 769080 fax: +44(0)1642 765848  
 e-mail: enquiries@implantsinternational.com  
 web: www.implantsinternational.com