



IMPLANTS

INTERNATIONAL

SPINAL & MAXILLO FACIAL



PEEK OPTIMA® **ANATOMIC**
Cervical Cage




PEEK OPTIMA[®] ANATOMIC Cervical Cage

The profile of our **ANATOMIC** cages facilitate intracanal navigation.

The **ANATOMIC** cages are 14.5 or 17mm wide and 13 or 14mm long respectively to ensure high stability and maximum bone end-plate interface.

The **ANATOMIC** cages are in 6 different heights from 5 mm to 10 mm.



Cervical Cages 14.5mm Wide

Ref.	Length	Height
SS7167	13 mm	5 mm
SS7168	13 mm	6 mm
SS7169	13 mm	7 mm
SS7170	13 mm	8 mm
SS7171	13 mm	9 mm
SS7172	13 mm	10 mm

Cervical Cages 17mm Wide

Ref.	Length	Height
SS7173	17 mm	5 mm
SS7174	17 mm	6 mm
SS7175	17 mm	7 mm
SS7176	17 mm	8 mm
SS7177	17 mm	9 mm
SS7178	17 mm	10 mm

March 10

DISTRIBUTED BY



Manufactured in the UK

Implants International Ltd
 71 Jay Avenue, Teesside Industrial Estate
 Thornaby-on-Tees, TS17 9LZ, United Kingdom
 tel: +44(0)1642 769080 fax: +44(0)1642 765848
 e-mail: enquiries@implantsinternational.com
 web: www.implantsinternational.com