



IMPLANTS

I N T E R N A T I O N A L

S P I N A L & M A X I L L O F A C I A L



PEEK OPTIMA[®] **ANATOMIC**
Intersomatic Lumbar Cage



PEEK OPTIMA® ANATOMIC Intersomatic Lumbar Cage



The anatomical ballistic profile of these cages facilitates intracanal navigation and minimises the risk of nerve damage.

The **ANATOMIC** cages are 11mm wide and 25mm long promoting high stability and maximum bone end-plate interface.

These **ANATOMIC** cages are in 7 different heights from 9mm to 15mm. They are available in 2 standard lordotic angles of 5° and 8°

X-ray markers in Tantalum

5° Lumbar Cages

Ref.	Lenght	Height
SS7154	25 mm	9 mm
SS7155	25 mm	10 mm
SS7156	25 mm	11 mm
SS7157	25 mm	12 mm
SS7158	25 mm	13 mm
SS7159	25 mm	14 mm
SS7160	25 mm	15 mm

8° Lumbar Cages

Ref.	Lenght	Height
SS7161	25 mm	10 mm
SS7162	25 mm	11 mm
SS7163	25 mm	12 mm
SS7164	25 mm	13 mm
SS7165	25 mm	14 mm
SS7166	25 mm	15 mm

March 10

DISTRIBUTED BY



Manufactured in the UK

Implants International Ltd
71 Jay Avenue, Teesside Industrial Estate
Thornaby-on-Tees, TS17 9LZ, United Kingdom
tel: +44(0)1642 769080 fax: +44(0)1642 765848
e-mail: enquiries@implantsinternational.com
web: www.implantsinternational.com