




IMPLANTS

INTERNATIONAL

TRAUMATOLOGY



COBRA™ CHS SYSTEM
TITANIUM



COBRA™ TITANIUM CHS SYSTEM

Cobra Plates

Titanium 6/4 ELI to ISO 5832 part 3

Type	Length	Cat #
Standard Plates		
2 Hole	61.2mm	TS 8050
3 Hole	75.0mm	TS 8051
4 Hole	88.2mm	TS 8052
6 Hole	115.2mm	TS 8053
8 Hole	142.2mm	TS 8054
10 Hole	169.2mm	TS 8055

Flanged Plates

Type	Length	Cat #
2 Hole	78.5mm	TS 8056
3 Hole	92.0mm	TS 8057
4 Hole	105.5mm	TS 8058
6 Hole	132.5mm	TS 8059
8 Hole	159.5mm	TS 8060
10 Hole	186.5mm	TS 8061

Cannulated Screw

Titanium 6/4 ELI to ISO 5832 part 3

Length	Ø	Cat #
30mm	6.5	TS 8086
35mm	6.5	TS 8087
40mm	6.5	TS 8088
45mm	6.5	TS 8089
50mm	6.5	TS 8090
55mm	6.5	TS 8091
60mm	6.5	TS 8092
65mm	6.5	TS 8093
70mm	6.5	TS 8094
75mm	6.5	TS 8095
80mm	6.5	TS 8096
85mm	6.5	TS 8097
90mm	6.5	TS 8098
95mm	6.5	TS 8099
100mm	6.5	TS 8100
105mm	6.5	TS 8101
110mm	6.5	TS 8102
115mm	6.5	TS 8103
120mm	6.5	TS 8104

Lag Screws & Barrels

Titanium 6/4 ELI to ISO 5832 part 3

Type	Length	Cat #
Short Barrel & Lag Screw		
135°	60 ES	TS 8062
135°	70 S	TS 8064
135°	80 M	TS 8066
135°	100 L	TS 8068
140°	60 ES	TS 8070
140°	70 S	TS 8072
140°	80 M	TS 8074
140°	100 L	TS 8076
145°	60 ES	TS 8078
145°	70 S	TS 8080
145°	80 M	TS 8082
145°	100 L	TS 8084

Standard Barrel & Lag Screw

135°	60 ES	TS 8063
135°	70 S	TS 8065
135°	80 M	TS 8067
135°	100 L	TS 8069
140°	60 ES	TS 8071
140°	70 S	TS 8073
140°	80 M	TS 8075
140°	100 L	TS 8077
145°	60 ES	TS 8079
145°	70 S	TS 8081
145°	80 M	TS 8083
145°	100 L	TS 8085

Washer

Titanium 6/4 ELI to ISO 5832 part 3

I.D.	O.D.	Cat #
7.0	12.5	TS 8105

Surgical Technique ST#032/05
2 Instrument Tray System

Diaphyseal Screw

Titanium 6/4 ELI to ISO 5832 part 3

Length	Ø	Cat #
16mm	4.5	TS 7772
18mm	4.5	TS 7773
20mm	4.5	TS 7774
22mm	4.5	TS 7775
24mm	4.5	TS 7776
25mm	4.5	TS 7742
26mm	4.5	TS 7777
28mm	4.5	TS 7778
30mm	4.5	TS 7743
32mm	4.5	TS 7779
34mm	4.5	TS 7780
35mm	4.5	TS 7744
36mm	4.5	TS 7781
38mm	4.5	TS 7782
40mm	4.5	TS 7745
42mm	4.5	TS 7783
44mm	4.5	TS 7784
45mm	4.5	TS 7746
46mm	4.5	TS 7785
48mm	4.5	TS 7786
50mm	4.5	TS 7747
52mm	4.5	TS 7787
54mm	4.5	TS 7788
56mm	4.5	TS 7789
58mm	4.5	TS 7790
60mm	4.5	TS 7748
65mm	4.5	TS 7791
70mm	4.5	TS 7792

March 10

DISTRIBUTED BY



0120



Manufactured in the UK

Implants International Ltd
71 Jay Avenue, Teesside Industrial Estate
Thornaby-on-Tees, TS17 9LZ, United Kingdom
tel: +44(0)1642 769080 fax: +44(0)1642 765848
e-mail: enquiries@implantsinternational.com
web: www.implantsinternational.com