




Flexible,
user friendly
Supra-Condylar
femoral nail





S-CON™ SUPRA-CONDYLAR NAIL SYSTEM

Designed as a flexible, user friendly Supra-Condylar femoral nail combining a number of unique features:

- A ballistic tip guides the nail into the medullary canal
- 6.5mm cancellous screws anchor the nail into the condylar block
- The guide finds both the proximal 5mm and the off-set 6.5mm distal holes
- The 12mm solid-section titanium alloy has superior strength
- Suitable for tibio-calcaneal fusions

S-CON™ Nail

Dia	Length	Cat
12mm	150mm	IM7425
12mm	200mm	IM7424
12mm	250mm	IM7464

Material: TIN Coated 6/4 Titanium

Diaphyseal Interlocking Screws

Dia	Length	Cat
5mm	25mm	TS7670
5mm	30mm	TS7502
5mm	35mm	TS7503
5mm	40mm	TS7504
5mm	45mm	TS7505
5mm	50mm	TS7506
5mm	60mm	TS7671

Material: TIN Coated 6/4 Titanium

Fully Threaded Cancellous Screws

Dia	Length	Cat
6.5mm	45mm	TS7524
6.5mm	50mm	TS7525
6.5mm	60mm	TS7526
6.5mm	65mm	TS7527
6.5mm	70mm	TS7528
6.5mm	80mm	TS7529
6.5mm	90mm	TS7529
6.5mm	100mm	TS7529

Material: TIN Coated 6/4 Titanium

Surgical Technique: Ref ST# 021/01
Instrument Kit: 1 Tray System

Sept 09



Manufactured in the UK

DISTRIBUTED BY

Implants International Ltd
71 Jay Avenue, Teesside Industrial Estate
Thornaby-on-Tees, TS17 9LZ, United Kingdom
tel: +44(0)1642 769080 fax: +44(0)1642 765848
e-mail: enquiries@implantsinternational.com
web: www.implantsinternational.com